

MED7[®]

GET BACK, FAST




TCR-150 *TECARTHERAPY*

Tecartherapy - short for "Transfer of Energy Capacitive and Resistive", is an advanced electrotherapy technique used for pain relief, tissue repair and muscle recovery. It involves the application of electromagnetic energy through capacitive and resistive modes, to stimulate tissues at various depths. Capacitive mode targets superficial tissues, while resistive mode targets deeper muscles and tendons.




MED7[®]

GET BACK, FAST



WHEN TO USE TCR-150?



**Degenerative soft tissue
diseases**

Tendon ruptures

Recent trauma


Tendinitis

Epicondylitis

Bursitis

MADE IN ITALY

TCR-150



HOW IS TCR-150 DIFFERENT?

- **COST EFFECTIVENESS:** TCR-150 is a feature packed device, which usually requires buying a much more expensive machine
- **FLEXIBILITY:** TCR-150 is highly portable and easy to transport when visiting clients or away fixtures/events
- **RELAIBILITY:** high quality Italian manufacture with an excellent reliability and durability track record

SPECIFICATIONS:

- New architecture software
- 7" capacitive touchscreen display
- 2 mm tempered glass screen
- Customizable programs
- Lightweight and portable
- Updates via USB ports
- Selectable protocols

MED7[®]
GET BACK, FAST



NEED MORE INFORMATION?



TECHNICAL DETAILS

- Tecnosix TCR150 allows you to work in 5 different modes
 - Manual Capacitive Resistive
 - Dynamic automatic mode
 - Bipolar mode
 - Preset Protocols
 - Custom Programs
- Emission frequencies: 0.5 MHz (+0 - 15% typical 448 kHz) - 1.0 MHz (+0 - 20% typical 0,8 MHz)
- Medical Power Supply: 300 W
- Maximum power: 150W
- Maximum absorbed power: 300 VA
- Power supply 90-240 Vac / 50-60 Hz
- Dimensions in mm: 170 x 315 x 390
- Weight: 5 kg c.a.

MED7 contact details are set out below. We welcome discussions about the appropriate solutions for your practice

We also have a range of funding solutions including rental programmes for most of our devices

www.med7.co.uk

info@med7.co.uk

(01342) 363 057

GET BACK FAST

16 Upper Woburn Place | London | WC1H 0AF